

# BIG PAW

Making wheelchair life safer





**BETTER SUPPORT**  
**BETTER PROTECTION**  
**BETTER CONTROL**

BIGPAW is designed to fit on top of the original foot plate, maintaining the wheelchairs integrity. No tools required for fitting.

BIGPAW delivers increased support, protection and control of your feet.

# THE ADVANTAGES

BIG  
PAW



Incorporating the rolled edges and the leg strap into the BIGPAW design, the following advantages are achieved:

**Support**

Maximum surface area, reducing pressure injuries.

**Protection**

NO castor- foot contact.

**Control**

Incorporated strap provides better rear leg control. Engage the split strap to prevent the legs from moving forward, caused by rough terrain or spasm.

# LEG STRAP



The leg strap ends are fitted with finger loops for easy opening. Both ends are stored behind the rear strap; when the terrain gets rough and extra leg control is required, simply bring both ends around in front of your legs where both ends will lock together via velcro.

# BIG PAW SPECS

## What do you get?

- BIGPAW constructed from 2 mm aluminum, anodized matt black.
- Attached rubber surface - suitable for bare feet, while providing excellent grip for shoes.
- 1 meter x 5cm strap (3 ft x 2 in)
- Mounting tape and alcohol wipes for surface preparation

## Sizes available

- 23 cm, 324 g (9 1/16 in, 11.4 oz)
- 24 cm, 335 g (9 29/64 in, 11.8 oz)
- 25 cm, 341 g (9 27/32 in, 12 oz)
- 27 cm, 352 g (10 5/8 in, 12.4 oz)
- 29 cm, 369 g (11 27/64 in, 13 oz)

All BIGPAW's are 15 cm (6 in) in depth.



Visit [www.bigpaw.com.au](http://www.bigpaw.com.au) for additional information, BIGPAW size template and attachment tutorials



[bluedog@git.com.au](mailto:bluedog@git.com.au)

bluedog  
designs